



Decontamination Procedures

Cleaning Solution

Hypochlorite solutionⁱ (one part Bleach to nine parts water). Bleach wipes, Lysol wipes or any other hypochlorite solution may be used.

Cleaning Procedure

Clean Twice Before Use and After Use!

Manikins & Equipment

Manikins and equipment (Mask, BVM, AED pads etc.) must be wiped down before student use or course begins and immediately after use or end of the course.

Tables

Tables must be wiped down at the end of training course.

ⁱ Per OSHA Safety Protocols